

**LISTA DI COMPATIBILITA' - COMPATIBILITY LIST**

Casa Costruttrice Brand	Modello - Model	Serie - Series										Connettore Connector
----------------------------	-----------------	----------------	--	--	--	--	--	--	--	--	--	-------------------------

Casa Costruttrice Brand	Modello - Model	ADULTO				PEDIATRICO		PEDIATRICO				Connettore Connector
		0060207				0060208		0060209				
<b>Philips Laerdal Medical</b>	Heartstream ForeRunner I AED model E, S, EM	●	●	●	●	-	-	-	-	-	-	
	Heartstart FR1 AED modelli-models 940010XX, 940020XX, 940030XX	●	●	●	●	-	-	-	-	-	-	
	Heartstream Heartstart FR2 AED modelli-models M3841A, M3861A	●	●	●	●	-	-	●	●	-	-	
	Heartstream HeartStart FR2+ AED, con modulo ECG-with ECG module, modelli-models M3840A, M3860A	●	●	●	●	-	-	●	●	-	-	
	CodeMaster modelli-models M1722A/B, M1723A/B, M1724A, M2475B (con cavo-with cable Philips M3507A e adattatore-and pad adapter 05-200)	●	●	●	●	●	●	-	-	-	●	
	HeartStart MRx modello-model M3535A (con cavo-with cable Philips M3508A)	●	●	●	●	●	●	-	-	-	●	
	Efficia DFM100	●	●	●	●	●	●	-	-	-	●	
<b>Agilent – Philips</b>	Heartstream/HeartStart XL modello-model M4735A (con cavo-with cable Philips M3508A)	●	●	●	●	●	●	-	-	-	●	
	Heartstream/HeartStart XLT modello-model M3500B (con cavo-with cable Philips M3508A)	●	●	●	●	●	●	-	-	-	●	
	Heartstream/HeartStart XL+ modello-model 861290 (con cavo-with cable Philips M3508A)	●	●	●	●	●	●	-	-	-	●	
	Heartstart 4000 modello-model M5500B (con cavo-with cable Philips M3508A)	●	●	●	●	●	●	-	-	-	●	

F7950PL, F7950WPL: riduzione di energia erogata al paziente: 150J sul defibrillatore diventano 50J sul paziente-reduced energy: defibrillator output of 150J is reduced to 50J



Casa Costruttrice Brand	Modello - Model	ADULTO				PEDIATRICO				Connettore Connector
		0060216				0060215				
<b>Zoll Medical Corp</b>	ZOLL PD 1200	●	●	●	●	●	●	●	●	
	ZOLL PD 1400	●	●	●	●	●	●	●	●	
	ZOLL PD 1600	●	●	●	●	●	●	●	●	
	ZOLL PD 1700	●	●	●	●	●	●	●	●	
	ZOLL PD 2000	●	●	●	●	●	●	●	●	
	ZOLL SERIE M-M SERIES	●	●	●	●	●	●	●	●	
	ZOLL SERIE E-E SERIES	●	●	●	●	●	●	●	●	
	ZOLL SERIE R-R SERIES	●	●	●	●	●	●	●	●	
ZOLL SERIE X-X SERIES	●	●	●	●	●	●	●	●		
<b>M&amp;B</b>	AED 7000	●	●	●	●	-	-	-	-	



		ADULTO 0060213				PEDIATRICO 0060214		PEDIATRICO 0060219 C/ Reductor				
<b>Physio-Control</b>	LIFEPAK 9 presa-socket QUICK-COMBO	●	●	●	●	●	●	-	-	●	●	
	LIFEPAK 10	●	●	●	●	●	●	-	-	●	●	
	LIFEPAK 12	●	●	●	●	●	●	-	-	●	●	
	LIFEPAK 15	●	●	●	●	●	●	-	-	●	●	
	LIFEPAK 20	●	●	●	●	●	●	-	-	●	●	
	LIFEPAK 500 AED con presa grigia-with grey socket	●	●	●	●	-	-	-	-	-	-	
	LIFEPAK 500 AED con presa rosa-with pink socket	●	●	●	●	-	-	●	●	-	-	
	LIFEPAK 1000 AED	●	●	●	●	-	-	●	●	-	-	
	LIFEPAK CR PLUS AED	-	●	-	●	-	-	●	●	-	-	
LIFEPAK EXPRESS AED	-	●	-	●	-	-	●	●	-	-		
<b>Bexen Cardio</b>	REANIBEX-200	●	●	(*)	(*)	-	-	-	-	-	-	
	REANIBEX-300	●	●	(*)	(*)	-	-	-	-	-	-	
	REANIBEX-500	●	●	(*)	(*)	-	-	-	-	-	-	
	REANIBEX-700	●	●	(*)	(*)	●	●	-	-	(*)	(*)	
	REANIBEX-800	●	●	(*)	(*)	●	●	-	-	(*)	(*)	
<b>Cardioline</b>	ELIFE 700	●	●	●	●	●	●	-	-	●	●	
<b>Mindray</b>	BENEHEART D3	●	●	●	●	●	●	-	-	●	●	
	BENEHEART D6	●	●	●	●	●	●	-	-	●	●	

F7952PL, F7952WPL: riduzione di energia erogata al paziente: 200J sul defibrillatore diventano 50J sul paziente-reduced energy: defibrillator output of 200J is reduced to 50J

(\*) Non per stimolazione - Not suitable for pace

		ADULTO 0060217				ADULTO 0060218			
<b>Agilent Hewlett Packard</b>	CODEMASTER 100	●	●	●	●	●	●	●	●
	CODEMASTER XE	●	●	●	●	●	●	●	●
	CODEMASTER XL	●	●	●	●	●	●	●	●
	CODEMASTER XL+	●	●	●	●	●	●	●	●
<b>Laerdal Medical</b>	HEARTSTREAM XL, XLT e-and XL+, con cavo-with cable M3507A	●	●	●	●	●	●	●	●
	HEARTSTART 4000	●	●	●	●	●	●	●	●

		ADULTO 0060200				PEDIATRICO 0060201		
<b>Schiller</b>	FRED	●	●	●	●	●	●	 
	FRED TELE	●	●	●	●	●	●	
	DEFIGARD DG 1002	●	●	●	●	●	●	
	DEFIGARD DG 3002	●	●	●	●	●	●	
<b>Esaote</b>	DEFIGARD DG 2000 EDOS	●	●	●	●	●	●	
	MDF SAED	●	●	●	●	●	●	
	MDF II	●	●	●	●	●	●	

		ADULTO 0060203				ADULTO 0060204		
<b>Nihon Kohden</b>	Cardiolife TEC-5500 series <b>con cavo-with cable JC-755V</b>	●	●	●	●	●	●	 
	Cardiolife TEC-5600 series <b>con cavo-with cable K342A (JC-855V)</b>	●	●	●	●	●	●	
	Cardiolife TEC-7600 series <b>con cavo-with cable JC-755V</b>	●	●	●	●	●	●	
	Cardiolife TEC-7700 series <b>con cavo-with cable JC-755V</b>	●	●	●	●	●	●	
	Cardiolife TEC-2202 <b>con cavo-with cable JC-755V</b>	●	●	●	●	●	●	
	Cardiolife TEC-2203 <b>con cavo-with cable JC-755V</b>	●	●	●	●	●	●	
	Cardiolife TEC-7531 <b>con cavo-with cable JC-755V</b>	●	●	●	●	●	●	
	Cardiolife TEC-8251 <b>con cavo-with cable JC-755V</b>	●	●	●	●	●	●	
	Cardiolife TEC-8300 series <b>con cavo-with cable K342A (JC-855V)</b>	●	●	●	●	●	●	
<b>Welch Allyn</b>	AED-10 (Adattatore pediatrico-paediatric adapter Welch Allyn 002173-U <sup>(2)</sup> )	●	●	●	●	●	●	
	AED-20 (Adattatore pediatrico-paediatric adapter Welch Allyn 002174-U <sup>(2)</sup> )	●	●	●	●	●	●	
	PIC-30, PIC-40, PIC-50	●	●	●	●	●	●	
<b>Artema Innomed</b>	CARDIOAID 360B	●	●	●	●	●	●	
	CARDIOAID 200B	●	●	●	●	●	●	
	CARDIOAID 200	●	●	●	●	●	●	
	CARDIOAID 100	●	●	●	●	●	●	
	CARDIOAID XL+	●	●	●	●	●	●	
	CARDIOAID XL	●	●	●	●	●	●	
<b>Artema S&amp;W</b>	CARDIO AID LS	●	●	●	●	●	●	
	CARDIO AID MC	●	●	●	●	●	●	
	CARDIO AID MC+	●	●	●	●	●	●	
<b>Dräger</b>	CARDIOLOG 2000	●	●	●	●	●	●	
<b>MRL</b>	LITE	●	●	●	●	●	●	
	LIFEQUEST	●	●	●	●	-	-	

		ADULTO 0060220				PEDIATRICO 0060221		
<b>Schiller</b>	FRED EASY	•	•	•	•			
	FRED EASY 2G	•	•	•	•			
	DEFIGARD DG 4000	•	•	•	•			
	DEFIGARD DG 5000	•	•	•	•			
	DEFIGARD DG 6002	•	•	•	•			
		ADULTO 0060210			PEDIATRICO 0060211	PEDIATRICO 0060212 C/ Reductor		
<b>CU Medical Systems, inc</b>	PARAMEDIC CU-ER1	•	•	-	•			
	PARAMEDIC CU-ER2, CU-ER3	•	•	•	-			
	PARAMEDIC CU-ER5 con cavo-with cable ER-5 OA03	•	•	•	-			
	LIFEGAIN CU-HD1	•	•	•	-			
		ADULTO 0060222			PEDIATRICO 0060223 C/ Reductor			
<b>CU Medical Systems, inc</b>	I – PAD NF1200	•	•	•	-	•	•	
<b>Cmos Drake</b>	FUTURA	•	•	•	•	-	-	
						ADULTO 0060202		
<b>Metrax Primedic</b>	HEARTSAVE AED (fino a-up to s/n 73944xxxxxx)					•	•	
	HEARTSAVE AED-M (fino a-up to s/n 73945xxxxxx)					•	•	
	HEARTSAVE 6 (fino a-up to s/n 73946xxxxxx)					•	•	
	HEARTSAVE 6S (fino a-up to s/n 73947xxxxxx)					•	•	
	DEFIMONITOR XD					•	•	

ADULTO 0060228			PEDIATRICO 0060229
-------------------	--	--	-----------------------

<b>Nihon Kohden</b>	Cardiolife TEC-2200 series <b>con cavo-with cable JC-765V</b>	●	●	●	●	
	Cardiolife TEC-2300 series <b>con cavo-with cable JC-765V</b>	●	●	●	●	
	Cardiolife TEC-5500 series <b>con cavo-with cable JC-765V</b>	●	●	●	●	
	Cardiolife TEC-5600 series <b>con cavo-with cable K342B (JC-865V)</b>	●	●	●	●	
	Cardiolife TEC-7600 series <b>con cavo-with cable JC-765V</b>	●	●	●	●	
	Cardiolife TEC-7700 series <b>con cavo-with cable JC-765V</b>	●	●	●	●	
	Cardiolife TEC-7531 <b>con cavo-with cable JC-765V</b>	●	●	●	●	
	Cardiolife TEC-8300 series <b>con cavo-with cable K342B (JC-865V)</b>	●	●	●	●	

ADULTO 0060205		PEDIATRICO 0060206
-------------------	--	-----------------------

<b>Marquette GE</b>	MARQUETTE 1250 FIRST RESPONDER	●	●	●	
	MARQUETTE 1500 FIRST RESPONDER	●	●	●	
	MARQUETTE RESPONDER 2000 (con cavo adattatore-with adapter cable)	●	●	●	
	MARQUETTE RESPONDER 3000 (con cavo adattatore-with adapter cable)	●	●	●	
<b>Esaote</b>	MDF	●	●	●	

ADULTO 0060226
-------------------

<b>Nihon Kohden</b>	TEC-8250	●	
	TEC-8251	●	

ADULTO 0060227	
-------------------	--

<b>Weinmann</b>	MEDUCORE EASY	●	-	
	MEDUCORE STANDARD	●	●	
<b>Cardiaid</b>	CARDIAID CT0207	●	-	

	ADULTO 0060230	PEDIATRICO 0060231	
--	-------------------	-----------------------	--

<b>Corpuls</b>	CORPULS 3		•	•		
----------------	-----------	--	---	---	--	--

Gli elettrodi non sono compatibili nel caso CORPULS 3 dotato del modular system Defibrillator-/Pacerunit **SLIM**-*Electrodes not compatible if CORPULS 3 is equipped with modular system Defibrillator-/Pacerunit SLIM*

ADULTO 0060232	PEDIATRICO 0060233
-------------------	-----------------------

<b>Defibtech</b>	DDU – 100 Lifeline AED	•	•		
	DDU – 120 Lifeline AUTO AED	•	•		

F7966WPL: riduzione di energia erogata al paziente: 150J sul defibrillatore diventano 50J sul paziente-*reduced energy: defibrillator output of 150J is reduced to 50J*

ADULTO 0060234	PEDIATRICO 0060235
-------------------	-----------------------

<b>Schiller</b>	FRED EasyPort	•	•	
	ARGUS PRO LifeCare	•	•	
	ARGUS PRO LifeCare 2	•	•	
	CARDIOVIT AT-101 Easy	•	•	

ADULTO 0060224	PEDIATRICO 0060225
-------------------	-----------------------

<b>Mediana</b>	HeartOn AED A10	•	•		
----------------	-----------------	---	---	--	--

ADULTO  
0060238


<b>GE</b>	RESPONDER 3000 (attacco diretto-direct connection)	●	●
	RESPONDER 2000 (attacco diretto-direct connection)	●	●
	RESPONDER 1100 (attacco diretto-direct connection)	●	●
	CARDIOSERV (attacco diretto-direct connection)	●	●
	DASH RESPONDER (attacco diretto-direct connection)	●	●



ADULTO  
0060239

PEDIATRICO  
0060240

<b>Metsis</b>	Life Point Pro AED	●	●
---------------	--------------------	---	---



ADULTO  
0060241

PEDIATRICO  
0060242

<b>Instramed</b>	INSTRAMED AED	●	●
------------------	---------------	---	---



ADULTO  
0060243

PEDIATRICO  
0060244

<b>Instramed</b>	INSTRAMED CARDIOMAX	●	●
------------------	---------------------	---	---



		ADULTO 0060245	PEDIATRICO 0060246	
Zoll Medical Corp	AED PLUS	•	•	
	AED PRO	•	•	

F7987W-F7987PW: Si prega di notare che manca il feedback CPR (profondità compressione & metronomo stimolazione) per queste piastre AED - *Please note there is no CPR feedback (compression depth & pace metronome) for these*