

ri-pen®

Diagnostic pupil penlight



Sustainable penlight for medical diagnosis with LED illumination

- ▶ Photobiologically validated for lighting of pupils according to EN 62471:2008
- ▶ UV free light through LED illumination
- ▶ Better vision due to a centered light beam
- ▶ Supplied singly or in packs of 6. (Assorted colours or six per color)
- ▶ Aluminium outer casing with disinfectable special lacquering
- ▶ Light aluminium casing with metal clip
- ▶ Simple on/off function through contact with the metal clip.
- ▶ Including 2 batteries type AAA.
- ▶ 2-year warranty

Technical data:

ri-pen®	3 V LED		
	5.300 Kelvin		
	40.000 Lux max.		
	Field of view approx. 40 mm diameter (distance 15 cm)		
Batteries:	2 x Alkaline AAA, 1,5 V		
No. 5075-526	ri-pen® black, pack of 1 pc.	No. 5075	ri-pen® black, pack of 6 pcs.
No. 5074-526	ri-pen® silver, pack of 1 pc.	No. 5074	ri-pen® silver, pack of 6 pcs.
No. 5077-526	ri-pen® red, pack of 1 pc.	No. 5077	ri-pen® red, pack of 6 pcs.
No. 5078-526	ri-pen® purple, pack of 1 pc.	No. 5078	ri-pen® purple, pack of 6 pcs.
No. 5071-526	ri-pen® blue, pack of 1 pc.	No. 5071	ri-pen® blue, pack of 6 pcs.
No. 5072-526	ri-pen® green, pack of 1 pc.	No. 5072	ri-pen® green, pack of 6 pcs.
		No. 5070	ri-pen® assorted colours, pack of 6 pcs.

Please visit riester.de for more information